

Kinectiv Technology Neck-Heads and Neck-Stems

Indicates "Combination is approved by Zimmer, Inc."

Straight Necks
(00-7848-002/044-00)
Anteverted Necks
(00-7848-002/044-01)
Retroverted Necks
(00-7848-002/044-01)

CoCr Femoral Heads, 12/14 Taper¹

DURASUL CoCr Femoral Head, +0 Length
METASUL Femoral Head +0 Length
METASUL LDH (Large Diameter Head)
METASUL LDH Head Adapter, +0 Length (12/14 taper)
VerSys CoCr Femoral Head, +0 Length

✓	✓	✓
✓	✓	✓
✓	✓	✓
✓	✓	✓
✓	✓	✓

Ceramic Femoral Heads, 12/14 Taper¹

BIOLOX-delta Femoral Head, +0 Length
BIOLOX-forte Femoral Head, +0 Length
BIOLOX-option Femoral Head, +0 Length
CERASUL Femoral Head System, +0 Length
SULOX Femoral Head, +0 Length ²
Alumina 12/14 Ceramic Heads, +0 Length

✓	✓	✓
✓	✓	✓
✓	✓	✓
✓	✓	✓
✓	✓ ²	✓ ²
✓	✓	✓

Stainless Steel Femoral Heads, 12/14 Taper¹

PROTASUL-S30 Femoral Head, +0 Length

✓	✓	✓
---	---	---

Bipolar Heads

Bipolar head (4404-22/26-XXX)
DURASUL Bi-polar head
Multipolar Bipolar Cup

✓	✓	✓
✓	✓	✓
✓	✓	✓

Femoral Stems^{1,2}

M/L Taper Femoral Stem with Kinectiv Technology
CLS Brevius Stem with Kinectiv Technology

✓	✓	✓
✓	✓	✓

1. Any Zimmer hip prosthesis with Kinectiv modular neck technology is designed to be used with a +0 head only. Use of longer heads could lead to fatigue failure of the stem or neck and use of shorter heads could lead to reduced range of motion and component impingement.

2. Due to component impingement, do not use SULOX ceramic heads with the following anteverted or retroverted modular necks:

- 00-7848-002-01 Modular Neck K1
- 00-7848-003-01 Modular Neck S1
- 00-7848-020-01 Modular Neck K2
- 00-7848-030-01 Modular Neck S2