

Articulating Combinations: Metal Femoral Heads and Polyethylene Articulation

Indicates "Combination is approved by Zimmer, Inc."
 Indicates "Combination is not approved by Zimmer, Inc."

Femoral Heads:	12/14 Femoral Head (VerSys system), except 22mm	12/14 Femoral Head (VerSys system), 22mm	6 Degree Femoral Head, except 22mm	6 Degree Femoral Head, 22mm	Continuum Bearing Head, Hi Carbon	Continuum Bearing Head, Metal Superpolished	Continuum J-Taper Head	CPT stainless steel head, except 22mm	CPT stainless steel head, 22mm	CoCr/DURASUL CoCr Head 28, 32, 36 and 40mm	CoCr/DURASUL CoCr Head 38 & 44mm	CoCr/DURASUL CoCr Head 22 & 26mm	CoCr/DURASUL CoCr Head 22.2mm	METASUL eccentric head, 28 & 32mm	METASUL standard head, 28, 32, 36 and 40mm	PROTASUL-S30 head, 22.2mm	PROTASUL-S30 head, 22mm	PROTASUL-S-30 head, 28 & 32mm	Testa Jupe head	Tete a Jupe Acier Inox Stainless Steel 316L Cone 12/14 D28	Tete a Jupe Inox HTA Cone 8/10 D28	Tete Chrome-Cobalt Cone 12/14 D28/48	TRIBOSUL head	Zimtron stainless steel femoral head, except 22mm	Zimtron stainless steel femoral head, 22mm	Monoblock Stems:	CDH stem PROTASUL-S30 (monoblock)	Moore monoblock stem	Straight stem PROTASUL-S30 (monoblock)	Thompson monoblock stem			
ALLOCLASSIC CSF DURASUL insert	✓	✓							✓																								
ALLOCLASSIC CSF insert	✓								✓				✓	✓	✓	✓	✓	✓	✓											✓	✓		
ALPHA DURASUL insert	✓	✓							✓																								
ALPHA PE insert	✓	✓							✓																								
APR PE insert	✓	✓							✓																								
BETA DURASUL insert	✓	✓							✓																								
BETA PE insert	✓	✓							✓																								
CLS SPOTORNO DURASUL insert*	✓	✓							✓																								
CLS SPOTORNO PE insert*	✓	✓							✓																								
CLW PE insert	✓	✓							✓																								
Continuum Hedrocel Acetabular Revision Shell Liner	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓																							
Continuum Hedrocel Replacement Cup Insert	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓																							
ESKA PE symmetric liner																																	
Epsilon DURASUL insert	✓	✓							✓																								
Epsilon DURASUL Constrained insert	✓	✓							✓																								
Epsilon PE Insert	✓	✓							✓																								
FITEK PE insert	✓	✓							✓																								
GAMMA DURASUL insert	✓	✓							✓																								
GAMMA PE insert	✓	✓							✓																								
HG II PE Liner, 22 & 26mm	✓	✓	✓						✓																								
HG II PE Liner, 28 & 32mm	✓	✓	✓						✓																								
HG II XLPE Liner, 22 & 26mm	✓	✓	✓						✓																								
HG II XLPE Liner, 28 & 32mm	✓	✓	✓						✓																								
Kappa PE insert	✓	✓							✓																								
LAMELLA PE insert	✓	✓							✓																								
Longevity IT Constrained Liner*	✓	✓							✓																								
Longevity IT Liners	✓	✓							✓																								
LOR PE insert	✓	✓							✓																								
St. Nabor PE insert	✓	✓							✓																								
STAFIT PE insert	✓	✓							✓																								
STANDARD PE insert	✓	✓							✓																								
TM Revision Cemented Constrained Liner*	✓	✓							✓																								
TM Revision Liner	✓	✓							✓																								
Trilogy Constrained Liner*	✓	✓							✓																								
Trilogy Longevity Constrained Liner*	✓	✓							✓																								
Trilogy Longevity Liner, 22 & 26mm	✓	✓							✓																								
Trilogy Longevity Liner, 28, 32, & 36mm	✓	✓							✓																								
Trilogy Longevity Liner, 40mm	✓	✓							✓																								
Trilogy PE Liner, 22 & 26mm	✓	✓							✓																								
Trilogy PE Liner, 28 & 32mm	✓	✓							✓																								
Zintra liner	✓	✓							✓																								

Continued on Page 2

* Legacy Zimmer skirted femoral heads are not approved by Zimmer for use with CLS SPOTORNO DURASUL inserts, CLS SPOTORNO PE inserts, Trilogy Constrained Liners, Trilogy Longevity Constrained Liners, TM Revision Cemented Constrained Liners or Longevity IT Constrained Liners.

This represents functional compatibility and does not necessarily indicate regulatory approval status for marketing of product pairs.

