

Articulating Combinations: Metal on Metal

- Indicates "Combination is approved by Zimmer, Inc."
- Indicates "Combination is not approved by Zimmer, Inc."

Metal femoral heads:

Continuum Bearing Head, Hi Carbon

DUROM Femoral component *

METASUL eccentric head, 28 & 32mm

Metasul LDH (Large Diameter Head)

METASUL standard head, 28 & 32mm

METASUL standard head, 36 & 40mm

Metal Inserts (Liners):

ALLOCLASSIC CSF METASUL insert

ALPHA METASUL insert

APR METASUL insert

BETA METASUL insert

CLS SPOTORNO METASUL insert

CLW METASUL insert

Epsilon METASUL insert

GAMMA METASUL insert

Kappa METASUL insert

LAMELLA METASUL insert

Metasul Taper Liner (8770)

St. Nabor METASUL insert

STANDARD METASUL insert

		✓		✓	
		✓		✓	
		✓		✓	
		✓		✓	
		✓		✓	
		✓		✓	
		✓		✓	
		✓		✓	
		✓		✓	
		✓		✓	
				✓	✓
		✓		✓	
		✓		✓	

Metal-lined Articulating Cups:

DUROM Acetabular Component *

Low Profile METASUL cup

Marburg METASUL cup

Morscher Press-fit METASUL cup

Weber METASUL cup

Zimmer MMC Cup *

	✓		✓		
		✓		✓	
		✓		✓	
		✓		✓	
		✓		✓	
	✓		✓		

*** NOTE: The Durom Femoral Component is not approved for use with the Durom Acetabular cup and Zimmer MMC Cup for total hip surface replacement Arthroplasty in the U.S.A.**