

Acetabular Reconstruction: Polyethylene Liners/Cups and Metal Cages/Rings

Indicates "Combination is approved by Zimmer, Inc."

Indicates "Combination is not approved by Zimmer, Inc."

	TM Revision Cage	Burch-Schneider Cage	Ganz Ring	Müller Ring	Wagner Reinforcement Cage	ZCA Cage	ZCA Roof Ring
Acora cup (cemented)	✓	✓	✓	✓			
Apollo All-poly (PE) Cemented cup	✓	✓	✓				
BETA DURASUL insert					✓		
BETA METASUL insert					✓		
BETA PE insert					✓		
Brunswick PE cup	✓	✓	✓				
Continuum Hedrocel Acetabular Revision Shell Liner							
Full Profile PE cup	✓	✓	✓				
HG II PE Liner							
HG II XLPE Liner							
Longevity IT Liner							
Low Profile DURASUL cup		✓	✓	✓			
Low Profile METASUL cup		✓	✓	✓			
Low Profile PE cup		✓	✓	✓			
TM Revision Cemented Constrained Liner	✓						
TM Revision Liner	✓						
TM Revision Shell (Liner/Cup/Cage Construct)	✓						
Trilogy AB Alumina Insert							
Trilogy Longevity Liner							
Trilogy PE Liner							
ZCA Cup						✓	✓