

### Glenosphere to Liner Combinations Reverse Shoulder System

- Indicates "Combination is approved by Zimmer Biomet"
- Indicates "Combination is not approved by Zimmer Biomet"

	TM Reverse Standard Poly Liner Size 36mm (+0mm, +3mm, +6mm )	TM Reverse 65° Retentive Poly Liner Size 36mm (+0mm, +3mm, +6mm)	Vivacit-E Liner Plus (0mm, 3mm, 6mm) Offset Size 36mm Diameter	Vivacit-E Retentive Liner Plus (0mm, 3mm, 6mm) Offset 65° Neck Angle Size 36mm	Anatomical Shoulder Inverse/Reverse Humeral PE-Inlay(0mm, +3mm, +6mm) Size 36mm	ARCOM XL 44-36 Standard Humeral Bearing (+0mm, +3mm)	ARCOM XL 44-36 Retentive + 3mm Humeral Bearing	E1 44-36 Standard Humeral Bearing (+0mm, +3mm)	E1 44-36 Retentive +3 Humeral Bearing	TM Reverse Standard Poly Liner Size 40mm (+0mm, +3mm, +6mm)	TM Reverse 65° Retentive Poly Liner Size 40 mm (+0mm, +3mm, +6mm)	Vivacit-E Liner Plus (0mm,+3mm, +6mm) Offset 40mm Diameter	Vivacit-E Retentive Liner Plus (0mm, 3mm, 6mm) Offset 65° Neck Angle Size 40mm	Anatomical Shoulder Inverse/Reverse Humeral PE-Inlay(0mm, +3mm, +6mm) Size 40mm	ArComXL Standard Humeral Bearing (44-41 mm) (+0mm, +3mm)	ArComXL Retentive +3 mm Humeral Bearing (44-41 mm)	E1 Standard Humeral Bearing (44-41 mm) (+0mm, +3mm)	E1 Retentive +3 mm Humeral Bearing (44-41 mm)
TM Reverse Glenosphere Size 36mm	✓	✓	✓	✓	✓													
Anatomical Shoulder Inverse/Reverse Glenoid Head Size 36mm	✓	✓	✓	✓	✓													
Versa-Dial Glenosphere Standard Size 36mm	✓	✓	✓	✓		✓	✓	✓	✓									
Versa-Dial Glenosphere +3mm Size 36mm	✓	✓	✓	✓		✓	✓	✓	✓									
Versa-Dial Glenosphere +6mm Size 36mm	✓	✓	✓	✓		✓	✓	✓	✓									
Versa-Dial Glenosphere Standard – Ti Size 36 mm						✓	✓	✓	✓									
Versa-Dial Glenosphere Standard – Ti +3 Size 36mm						✓	✓	✓	✓									
Versa-Dial Glenosphere Standard – Ti +6 Size 36mm						✓	✓	✓	✓									
TM Reverse Glenosphere Size 40mm										✓	✓	✓	✓	✓				
Anatomical Shoulder Inverse/Reverse Glenoid Head Size 40mm										✓	✓	✓	✓	✓				
Versa-Dial Glenosphere Standard Size 41mm															✓	✓	✓	✓
Versa-Dial Glenosphere +3mm Size 41mm															✓	✓	✓	✓
Versa-Dial Glenosphere +6mm Size 41mm															✓	✓	✓	✓
Versa-Dial Glenosphere Standard – Ti Size 41 mm															✓	✓	✓	✓
Versa-Dial Glenosphere Standard – Ti +3 Size 41mm															✓	✓	✓	✓
Versa-Dial Glenosphere Standard – Ti +6 Size 41mm															✓	✓	✓	✓